

#### H301-ALPHA3-LARGE

H301-ALPHA3-LARGE is compatible with PhaseView Alpha<sup>3</sup> Light Sheet microscope.

IST 2068\_REV01

# H301-ALPHA3-LARGE



#### Index

1.	Components and dimensions	. 3
2.	Installation of the Chamber around the sample holder	. 4
3.	Connection of the Gas Supply	. 6
4.	Working with perfusion	. 6
5.	Insertion of the Sample Feedback Temperature Sensor	. 7
6.	Cleaning	. 8



# 1. Components and dimensions

H301-ALPHA3-LARGE includes the following components:

- Objective Collar, compatible with water immersion objective 21 mm diameter.
- Sample Arm Cover
- Chamber Lid
- Chamber Base




Figure 1. H301-ALPHA3-LARGE - Components and dimensions.


# H301-ALPHA3-LARGE



### 2. Installation of the Chamber around the sample holder.

To install the H301-ALPHA3-LARGE around the sample holder, please follow the steps below:

- 1. Insert the chamber base by sliding into place as shown in image 1 and 2 in Figure 2. **Note:** Please pay attention to the sample holder glass.
- 2. Tighten the captive screw M2 to fix the chamber base around the sample holder by using the provided 1.5 mm hexagonal key (see image 2 in Figure 2).
- 3. Place the magnetic lid on the chamber base (see image 3 in Figure 2).
- 4. Insert the sample arm cover (see image 4 in Figure 2), and close it (see image 5 in Figure 2). It allows to reduce the brightness into the sample.
- 5. Place the magnetic objective heater in place as shown in image 6 in Figure 2.





#### 3. Connection of the Gas Supply

A single silicon tubing carries output gas from the Okolab Gas Controller to H301-ALPHA3-LARGE. Silicon tubing connects to a gas input - brass opening - located on a corner of the H301-ALPHA3-LARGE (see Figure 3). Connect by gently pushing silicon tubing onto brass opening.



Figure 3. Connection with gas supply.

## 4. Working with perfusion

The Chamber H301-ALPHA3-LARGE features 2 perfusion holes on the sides for the insertion of perfusion tubing up to 1 mm in outer diameter. Follow the instruction below to work with perfusion:

- 6. Use the supplied Stainless Steel Perfusion Tubes. Attach a Silicon tubes O.D. 1,6mm I.D. 0,8mm (not supplied) to the Stainless Steel Perfusion Tubes (see image 1 in Figure 4).
- 7. Insert the Stainless Steel Perfusion Tubes in the perfusion holes (see images 2 and 3 in Figure 4).

  Insert a Silicon tubes O.D. 1,6mm I.D. 0,8mm (not supplied) to the free end of the Stainless Steel Perfusion Tubes (see image 4 and 5 in Figure 4).




Figure 4. Working with Perfusion

#### 5. Insertion of the Sample Feedback Temperature Sensor

Insert the Sample Feedback Temperature Sensor through the dedicated opening located in the H301-ALPHA3-LARGE (see Figure 5).

To read the temperature of the culture medium and realize a self-calibration of the chamber, fix the Temperature Sensor tip with adhesive tape on the sample holder inside the H301-ALPHA3-LARGE. For more information about the self-calibration, please refer to the User Manual of the Okolab Controller in use.






Figure 5. Insertion of the temperature sensor inside the chamber.

### 6. Cleaning

- Unplug all electrical cables.
- Remove the chamber lid from the chamber main body, and keep it separate from the chamber base while the chamber cools down.
- To clean the body of the chamber, wipe with a soft micro-fiber cloth. For stubborn smudges, you can damp
  the soft micro-fiber cloth with ethyl alcohol (product code UN1170). Do not put any liquid directly on the
  chamber.

H301-ALPHA3-LARGE\_User Manual Rev 01\_Published in April 2021