



# Okolab Sample Holders & Stage Adapters: Features and Benefits



## Cost Effective



- Okolab Universal Sample Inserts and Stage Adapters are affordably priced.
- No need to double up on fixed and live cell specimen holders.

## Quality Built



- Our Inserts are robustly built to withstand daily use in and out of our incubators.
- They are easy to use and designed with ergonomics in mind.

## Versatile



- Get more value from your Okolab Inserts by using them for all on-stage applications, even without an incubation chamber - simply add a multiwell plate adapter that fits your stage.
- Future-proof your setup: in case you add incubation later on and already have our Universal Inserts, they're compatible with nearly all Okolab chamber models.
- Use our four-slide holder 4xGS-INSERT to transform your inverted system to a slide scanner!



### 4xGS-INSERT

All INSERTS fit directly inside a multiwell plate holder. Use your own, or use one from Okolab, shown below.

### K-FRAME-MW-HOLDER

Hosts a multiwell plate, or any Okolab INSERT.

*Fits most K-FRAME stages, including:  
ASI/Marzhauer/Zaber.*



 [usa@oko-lab.com](mailto:usa@oko-lab.com)

 <https://oko-lab.com/>



# Available Universal Inserts & Adapters:



GS35-INSERT



1xLABTEK-INSERT



1x60-INSERT



2xLABTEK-II-INSERT

1x35-INSERT

4x35-INSERT

1xLABTEK-II-INSERT

1xGS-INSERT

2x35-INSERT

LABTEK-35-INSERT

2xGS-INSERT

LABTEK-II-35-INSERT

2xLABTEK-INSERT

LABTEK-60-INSERT

2x60-INSERT

LABTEK-II-60-INSERT

3xGS-INSERT

GS60-INSERT


4xGS-INSERT

6035-INSERT

*Our wide variety of INSERTS accommodate nearly all commercial culture dishes:*

*GS = 75x25mm | LABTEK = 57x25mm | LABTEK-II = 55x24mm*

**NEW!**



## ADJUSTABLE-INSERT

Perfect when an adjustable-style holder is preferred for imaging FIXED samples, even when the Okolab chamber is kept installed on the stage. ADJUSTABLE-INSERT uses two smoothly running brackets with a variable clamping range. It can host round dishes with diameter from 35 to 60mm, and slides with dimensions 25-54 x 24-92mm. Fits inside Okolab H203/H303/UNO chambers, or in any multiwell plate holder.

### Adapters

#### K-FRAME-MW-HOLDER:

Holds a multiwell plate or Okolab INSERTS in stages with K-Frame opening.

Adapters to fit other stage models are available.

Contact us for a compatibility check or with any questions.

#### Contact Our Sales Team Today!

Lara Petrak

Director of Sales, North America

Tyler Machovina

North American Sales Manager



petrak@oko-lab.com

machovina@oko-lab.com



650-410-0756

650-228-6154